

ASAHI CHIKAI Family



ACT ONE: The multi-wire winding offers 1:1 rotational transmission, shape memory and durability.



ASAHI CHIKAI 008

Micro guide wire compatible with small lumen micro catheters

Precision cases – super selective - distal treatments



ASAHI CHIKAI X 010

* Newest guide wire technology without ACT ONE

Very soft, distal shaft with good rotational stability even in distal vessels

Precision cases – super selective - distal treatments



ASAHI CHIKAI 10

Micro guide wire with high support power and stability

Precision cases - treatment of small vessels



ASAHI CHIKAI

Standard working wire provides support and stability

Standard cases - low complexity



ASAHI CHIKAI black

Polyurethane coating for higher gliding performance in complex anatomy

Tortuous vessels - high complexity



ASAHI CHIKAI black 14 Soft Tip

43% softer tip than CHIKAI black

Atraumatic tip - high complexity



ASAHI CHIKAI black 18

Tapered shape (distal 0.45mm (0.018 inches)-proximal 0.36 mm (0.014 inches)). Especially suitable for larger lumen catheters

Large vessels - high supportability



ASAHI CHIKAI black 18 Soft Tip

30% softer tip than CHIKAI 18

Atraumatic tip - Large vessel

