





ASAHI CHIKAI FAMILY

PRODUCT LINE UP



YOUR DREAMS. WOVEN TOGETHER.

ASAHI CHIKAI Family



Characteristics

1:1 torque response

ACT ONE technology provides a 1:1 torque response, in tortuous anatomies.

Shape of durable distal tip

The shape of the distal end is durable, especially for the following shapes 90° angled and "round curve".



※ Available for CHIKAI black



※ Available for CHIKAI, CHIKAI black, CHIKAI black 18, CHIKAI black 14/18 soft tip

ASAHI CHIKAI black Family

Neurovascular Guide Wire

Providing excellent support and less friction due to polymer jacket and hydrophilic coating for smooth delivery of microcatheter



ASAHI CHIKAI black (200cm)

0.36mm (0.014") guide wire with enhanced lubricity

0.014"




- Tip radiopacity 3cm
- Polymer jacket length 150cm
- SLIP-COAT® hydrophilic coating 170cm
- Complex cases and targeting distal anatomy
- Ease distal aspiration catheter placement

0.36mm (0.014")

ASAHI CHIKAI black 14 Soft Tip (200cm)

Safe manipulation

0.014"



- Tip radiopacity 3cm
- Polymer jacket length 150cm
- SLIP-COAT® hydrophilic coating 170cm
- Complex cases and targeting distal anatomy
- Ease distal aspiration catheter placement
- 43% softer tip than CHIKAI black
- Softer tip enhances safety in challenging anatomies

0.36mm (0.014")

ASAHI CHIKAI black 18 (200cm)

Speciality guide wire

0.018"



- Tip radiopacity 5cm
- Polymer jacket length 150cm
- SLIP-COAT® hydrophilic coating 170cm
- Enough support to deliver large ID microcatheter
- Reverse tapered shaft design to avoid ledge effect
- Ease clot retriever delivery and positioning

0.41mm (0.016")

0.36mm (0.014")

ASAHI CHIKAI black 18 Soft Tip (200cm)

Speciality guide wire

0.018"



- Tip radiopacity 5cm
- Polymer jacket length 150cm
- SLIP-COAT® hydrophilic coating 170cm
- Enough support to deliver large ID microcatheter
- Reverse tapered shaft design to avoid ledge effect
- Ease clot retriever delivery and positioning

0.41mm (0.016")


0.36mm (0.014")

- 30% softer tip than CHIKAI black 18

ASAHI CHIKAI Family


Neurovascular Guide Wire

ASAHI CHIKAI family are designed for an optimal balance between tip flexibility and microcatheter support



ASAHI CHIKAI (200cm, 300cm)


0.014" workhorse guide wire



- Tip radiopacity 5cm
- SLIP-COAT® hydrophilic coating 170cm
- Well-balanced workhorse wire

ASAHI CHIKAI 008 (200cm)


Soft flexible tip & distal body



- Tip radiopacity 9cm
- SLIP-COAT® hydrophilic coating 180cm
- Distal anatomy / small-sized vessels
- Compatible with microcatheters with small inner lumen

ASAHI CHIKAI 10 (200cm, 300cm)

Balanced tip and higher body support



- Tip radiopacity 3cm
- SLIP-COAT® hydrophilic coating 170cm
- Mid to distal anatomy / larger-sized vessels
- Compatible with microcatheters with relatively larger inner lumen


ASAHI CHIKAI X 010

Neurovascular Guide Wire



ASAHI CHIKAI X 010 (200cm)

Superior torque response in distal vessels




- Tip radiopacity 16cm
- Polymer jacket length 150cm
- SLIP-COAT® hydrophilic coating 180cm®

New Tip Core Design

Square Core-easy shaping, long lasting shape retention



Square Core
(ASAHI CHIKAI X 010)



Round Core
(Other ASAHI CHIKAI wires)



Structure & Order Information

Product Name	Product Code	Outer Diameter	Coating	Usable Length	Coil Length	Radiopaque Length	Tip Shape
ASAHI CHIKAI 008	WAIN-CKI-008-200	0.20/0.26mm (0.008/0.010inch)	Hydrophilic Coating (SLIP-COAT®) 180cm	200cm	9cm	9cm	Straight
ASAHI CHIKAI 10	WAIN-CKI-10-200	0.26mm (0.010inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	9.5cm	3cm	Straight
	WAIN-CKI-10-300	0.26mm (0.010inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	300cm	9.5cm	3cm	Straight
ASAHI CHIKAI X 010	AIN-CKX-10-200-R	0.26mm (0.010inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 180cm	200cm	16cm	16cm	Re-shapeable round curve (25° angle)
ASAHI CHIKAI	WAIN-CKI-200	0.36mm (0.014inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	30cm	5cm	Straight
	WAIN-CKI-200-RC	0.36mm (0.014inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	30m	5cm	Round curve 
	WAIN-CKI-300	0.36mm (0.014inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	300cm	30cm	5m	Straight
	WAIN-CKI-300-RC	0.36mm (0.014inch)	Hydrophilic Coating (SLIP-COAT®) 170cm	300cm	30cm	5cm	Round curve 
ASAHI CHIKAI black	WAIN-CKI-200-BS	0.36mm (0.014inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	30cm	3cm	Round curve 
	WAIN-CKI-200-BA	0.36mm (0.014inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	30cm	3cm	Angled 90° 
ASAHI CHIKAI black 14 soft tip	AIN-CKI-200-B-SFT	0.36mm (0.014inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	30cm	3cm	Round curve 
ASAHI CHIKAI black 18	WAIN-CKI-18-200-BS	0.45/0.36mm (0.018/0.014inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	34cm	5cm	Round curve 
ASAHI CHIKAI black 18 soft tip	AIN-CKI-18-200-SFT	0.45/0.36mm (0.018/0.014inch)	Polymer Jacket 150cm / Hydrophilic Coating (SLIP-COAT®) 170cm	200cm	34cm	5cm	Round curve 

Manufactured by:

 ASAHI INTECC

ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan

Contact phone number: +81-561-86-9101

www.asahi-intecc.com

Contact:

 ASAHI INTECC EUROPE B.V.

ASAHI INTECC EUROPE B.V.

Strawinskylaan 1461 (WTC Tower Ten - 14th floor)

1077 XX Amsterdam, The Netherlands

 0344



for Android



for iOS

ASAHI product information

K23234_B_LC_2.0