




ASAHI FUBUKI Guide Catheter


PRODUCT LINE UP



YOUR DREAMS. WOVEN TOGETHER.

ASAHI FUBUKI

A neurovascular guide catheter in a wide range of sizes.




Features

ASAHI FUBUKI Guide Catheter 6Fr, 7Fr, 8Fr


- **Flexibility**
 - Soft tip and flexible distal shaft minimizes vessel injury.
- **Pushability**
 - Balanced shaft with increased rigidity provides stable deployment and optimal pushability.
- **Lubricity**
 - Distal hydrophilic coating allows smooth navigation through aortic arch and vessel ostium.
- **Visibility**
 - Radiopaque marker tip & shaft enhance visibility provided by a platinum marker on the tip and radiopaque material in the shaft.

ASAHI FUBUKI

Fr	I.D.	O.D.
6	1.80 mm (0.071")	2.09 mm (6Fr)
7	2.05 mm (0.081")	2.40 mm (7Fr)
8	2.28 mm (0.090")	2.70 mm (8Fr)



Hydrophilic Coating	
6Fr	15cm
7Fr	15cm
8Fr	5cm






ASAHI FUBUKI Guide Catheter

ASAHI FUBUKI Guide Catheter 6Fr, 7Fr, 8Fr

- Balanced shaft rigidity for traceability and support
- Large inner lumen
- Good visibility

Use for various indications:




Structure & Ordering Information

Product Name	Product Code	Usable Length	Coating Length	Tip Shape
ASAHI FUBUKI 6Fr	WAIN-FBK-6A	90 cm	15 cm	Angle
	WAIN-FBK-6S	90 cm	15 cm	Straight
	WAIN-FBK-6A80	80 cm	15 cm	Angle
	WAIN-FBK-6S80	80 cm	15 cm	Straight
	WAIN-FBK-6AL	100 cm	15 cm	Angle
	WAIN-FBK-6SL	100 cm	15 cm	Straight
	WAIN-FBK-6A110	110 cm	15 cm	Angle
	WAIN-FBK-6S110	110 cm	15 cm	Straight
ASAHI FUBUKI 7Fr	WAIN-FBK-7A	90 cm	15 cm	Angle
	WAIN-FBK-7S	90 cm	15 cm	Straight
	WAIN-FBK-7A80	80 cm	15 cm	Angle
	WAIN-FBK-7S80	80 cm	15 cm	Straight
	WAIN-FBK-7AL	100 cm	15 cm	Angle
	WAIN-FBK-7SL	100 cm	15 cm	Straight
	WAIN-FBK-7A110	110 cm	15 cm	Angle
	WAIN-FBK-7S110	110 cm	15 cm	Straight
ASAHI FUBUKI 8Fr	WAIN-FBK-8S	90 cm	5 cm	Straight
	WAIN-FBK-8S80	80 cm	5 cm	Straight
	WAIN-FBK-8SL	100 cm	5 cm	Straight
	WAIN-FBK-8S110	110 cm	5 cm	Straight

ASAHI FUBUKI - Straight 

ASAHI FUBUKI - 45° Angle 



© 2021 ASAHI INTECC CO., LTD.
 "ASAHI" and "FUBUKI" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.
 Some materials / graphics were provided by Vascular Perspectives - www.vascularperspectives.com.

Manufactured by :

 **ASAHI INTECC**

ASAHI INTECC CO.,LTD.

Strawinskylaan 1461
 (WTC Tower Ten - 14th floor)
 1077 XX

<http://www.asahi-intecc.eu>

P21064_B_CE+_1.1