



PCI Guide Wires



Frontline		CLIP COLORS	TIP LOAD	
	ASAHI SION blue Well balanced workhorse		0.5gf	
	ASAHI SION Tracks tortuous vessels		0.7gf	
	ASAHI SION black Sliding performance		0.8gf	
Support				
	ASAHI SION blue ES Support		0.5gf	
	Grand Slam Support		0.7gf	
Complex & CTO				
	Fielder XT-R Navigates small vessels and tight lesions		0.6gf	
	Fielder XT-A Penetrates tapered proximal caps		1gf	
	ASAHI Miracle Neo 3 Controlled intraluminal crossing		3gf	
	ASAHI Gladius EX Controlled sliding		3gf	
	ASAHI Gladius MG Narrow loop technique		3gf	
	ASAHI Gaia Next Family Torque and tip flexibility		2gf 4gf 6gf	
	ASAHI Conianza Pro / Conquest Pro Penetration technique		9gf 12gf	
Retrograde				
	ASAHI SUOH 03 Crossability in angulated and tortuous collaterals		0.3gf	
	ASAHI RG3 Externalization		3gf	

Coatings Key: *Hydrophilic Hydrophobic (Silicone) Polymer Jacket & Hydrophilic Coating

©2025 ASAHI INTECC CO., LTD. "ASAHI", "SION", "Grand Slam", "Fielder", and "ASAHI Miracle" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries. Some materials / graphics were created in partnership with IP MED Group - www.ipmed.com
The tip load values of our guide wires are the values assuming the force that the tip of the guide wire pushes the lesion on the basis of company standardized test. There are individual differences in values of actual devices.
K25210_P_10_1.0