





# MSK EMBOLISATION PORTFOLIO




**Osteoarthritis**




**Sports Elbow**



**Frozen Shoulder**




**Wrist Pain**



## ASAHI Veloute Embolisation Microcatheter

- Low Profile and high trackability
- Designed for ultra-selective embolisation in tortuous anatomy



## ASAHI Meister 16

0.016" Guidewire

- High torque control for vessel selectivity
- Polymer jacket provides lubricity and trackability



# MSK EMBOLISATION PORTFOLIO

## Technical specs & ordering information

Product	Catalog No.	Usable Length	MC Length	O.D. of Distal Shaft/ O.D. of Proximal Shaft	I.D. of Distal Shaft/ I.D. of Proximal Shaft	Max. Guide Wire	Hydrophilic Coating Length	Tip Shape	Max. Pressure
ASAHI Veloute	VEL105-16Y	105 cm	112 cm	0.58/0.94 mm (1.7/2.8 Fr)	0.43/0.59 mm (0.017/0.023 inch)	0.41 mm (0.016 inch)	65 cm	Straight	6895 kPa (1000 psi)
	VEL105-16AY	105 cm	112 cm				65 cm	Angled	
	VEL125-16Y	125 cm	132 cm				65 cm	Straight	
	VEL125-16AY	125 cm	132 cm				65 cm	Angled	
	VEL150-16Y	150 cm	157 cm				100 cm	Straight	
	VEL150-16AY	150 cm	157 cm				100 cm	Angled	


Product	Catalog No.	Usable Length	Shaft O.D.	Radiopaque Length	Polymer Jacket Length	Hydrophilic Coating Length	Tip Shape
ASAHI Meister 16	WAMS-165-1645	165 cm	0.41 mm (0.016 inch)	5 cm	160 cm	160 cm	Angled 45°
	WAMS-180-1645	180 cm	0.41 mm (0.016 inch)	5 cm	175 cm	175 cm	Angled 45°
	WAMS-165-16WA	165 cm	0.41 mm (0.016 inch)	5 cm	160 cm	160 cm	Angled Double
	WAMS-180-16WA	180 cm	0.41 mm (0.016 inch)	5 cm	175 cm	175 cm	Angled Double
	WAMS-165-16ST	165 cm	0.41 mm (0.016 inch)	5 cm	160 cm	160 cm	Round Curve
	WAMS-180-16ST	180 cm	0.41 mm (0.016 inch)	5 cm	175 cm	175 cm	Round Curve

### Accessories



VEL105-16Y / VEL105-16AY  
VEL125-16Y / VEL125-16AY  
VEL150-16Y / VEL150-16AY

1. Stylet
2. 1cc Syringe
3. Y Connector



Angled 45°



Angled double



Round curve



ASAHI INTECC Europe B.V.

WTC Tower Ten, 14th floor, Strawinskyalaan 1461, 1077 XX Amsterdam, The Netherlands  
Contact phone number: +31-20-794-0642

K25229 P LC 1.0

©2026 ASAHI INTECC CO., LTD.

"ASAHI", "Veloute", and "ASAHI Meister" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries. All other trademarks are the property of their respective owners. Coated with SLIP-COAT® coating.

Some materials / graphics were created in partnership with VP MED Group - www.vpmed.com.